



	Instructions for use For Antiseptic Tulle Gras Dressing	Doc. No.: IU/24
		Rev. No.: 01
		Rev. Date: 13SEP2024

Descriptions:

The Antiseptic Tulle Gras Dressing is an antimicrobial dressing and made from an open weave gauze impregnated with White Soft Paraffin containing 0.5% w/w Chlorhexidine Acetate B.P. Dressing facilitate the passage of viscous exudate into a secondary absorbent dressing. It is presented sterile and interleaved for easy aseptic application. Chlorhexidine Acetate is a proven antiseptic which is active against a wide range of Gram-positive and Gram-negative organisms. It is indicated for use as a low adherent wound contact layer- on wound which are or maybe infected. For already infected wounds it can be used in conjunction with systemic antibacterial agents. It also allows easy drainage of the wound.

1. Indications:

Antiseptic Tulle Gras Dressing is an additive, topical, medicated dressing for:

- Minor open wounds.
- Minor burns and scalds.
- Laceration, abrasion and other skin loss wound.

2. Contraindications:

These dressings should not be used in known cases of sensitivity or allergy to Chlorhexidine Acetate.

3. Precautions:

- It should not be used on more than 10% of the body area.
- Dermatitis, contact sensitivity or photosensitivity may rarely occur.
- Avoid contact with eyes, the middle ear, meninges and brain.
- It is incompatible with soaps, other anionic materials, as well as potassium iodide.
- The use of Antiseptic Tulle Gras Dressing with oxidizing wound cleaning agents should be avoided.

4. Presentation

Each dressing (Including the roll and box) is wrapped with vegetable parchment paper in an individual foil pouch or box, sterilized by gamma radiation.

5. Packing and Labeling:

The Dressings are supplied in sterile condition; packed in individual aluminum foil pouches, polypropylene boxes and kept inside suitable size of outer box. Sterilization is carried out using the gamma irradiation process indicated by **STERILE R** on the label.


6. Storage condition:

Store flat in a dry place at a temperature not exceeding 25°C.
Protect devices from exposure to direct sunlight, heat and rain.



	Instructions for use For Antiseptic Tulle Gras Dressing	Doc. No.: IU/24
		Rev. No.: 01
		Rev. Date: 13SEP2024

7. Symbols Used in IFUs / Labels and Packaging Materials

Symbol	Definition	Symbol	Definition	Symbol	Definition
	Single use (Do not re-use)		Batch Number		Keep away from heat /sunlight and radioactive sources
	Date of Manufacture YYYY-MM		Manufactured by		Use by Date (Date of Expiry) YYYY-MM
	Caution: check for specific warnings or precautions		Do not re-sterilize		Do not use if package is opened or damaged
	Sterilized by Gamma Sterilization process		Consult instructions for use		Avoid moisture or water contact
	Recycle		Code Number / Part No.		To be sold only against prescription

8. Further information:

For further information concerning use of these devices, please check with Adler Customer Service at the address given herein or e-mail to adler-customer.care@adler-healthcare.com.

Manufactured by:



Adler Healthcare Pvt. Ltd.

Plot No. A-1, MIDC Sadavali, Tal. Sangameshwar, Dist. Ratnagiri, Maharashtra State, INDIA,
 PIN 415804. Ph: +91 (0)8055136000, 8055146000