




**Femoral Fracture Nail
Retrograde - Ti**



ATLAS Ti Femoral Fracture Nail Retrograde - Specifications

Material	Ti6AL4V
Diameters	10.0mm & 11.5mm
Lengths	28cm - 44cm
Nail Colour	Blue
Distal Diameter Driving End	12.0mm
Proximal Diameter Non-Driving End	10.0mm & 11.5mm
Locking Bolt Threads	M8 x 1
Screw Diameter	5.0mm
Screw Colour	Gold
Screw Lengths	20mm - 110mm
Screw Socket	Hex 4.75mm
Sterilization	Radiation
Distal Locking Driving End	
Static Lock Locations/Orientations	15mm/ML Can be Locked with TFN/FFN Nail Cap 30mm/25° Threaded with bushing 40mm/25° Threaded with bushing
Proximal Locking Non-Driving End	
Static Lock Locations/Orientations	15mm/AP 35mm/AP
AP Bow Radius	2000mm
AP Bow Location	Starts 100mm from Driving End



ATLAS Ti FFN Retrograde, Sterile

Code No.		Length cm
Ø 10mm	Ø 11mm	
I0146.1028s	I0147.1128s	28
I0146.1030s	I0147.1130s	30
I0146.1032s	I0147.1132s	32
I0146.1034s	I0147.1134s	34
I0146.1036s	I0147.1136s	36
I0146.1038s	I0147.1138s	38
I0146.1040s	I0147.1140s	40
I0146.1042*	I0147.1142*	42
I0146.1044*	I0147.1144*	44

*42 cm & 44 cm are supplied Non-Sterile

ATLAS Ti TFN/FFN Closing Cap

Code No.	Length mm
I0303.00	0
I0303.05	5
I0303.10	10

ATLAS Ti Locking Screw Ø5.0mm

Code No.	Length mm
Ø 5.0mm Color: Gold	
I0211.5020	20.0
I0211.5023	22.5
I0211.5025	25.0
I0211.5028	27.5
I0211.5030	30.0
I0211.5033	32.5
I0211.5035	35.0
I0211.5038	37.5
I0211.5040	40.0
I0211.5043	42.5
I0211.5045	45.0
I0211.5048	47.5
I0211.5050	50.0
I0211.5053	52.5
I0211.5055	55.0
I0211.5058	57.5
I0211.5060	60.0
I0211.5063	62.5
I0211.5065	65.0
I0211.5068	67.5
I0211.5070	70.0
I0211.5073	72.5
I0211.5075	75.0
I0211.5078	77.5
I0211.5080	80.0
I0211.5085	85.0
I0211.5090	90.0
I0211.5095	95.0
I0211.5100	100.0
I0211.5105	105.0
I0211.5110	110.0

Illustrations not to scale. Specifications subject to change without notice.



® Manufactured & marketed by

ADLER HEALTHCARE PVT. LTD.

Plot # A-1, MIDC, Sadavali (Devrukh) Tal. Sangameshwar, Dist. Ratnagiri. Pin - 415804. India.

Tel: 8055136000 / 8055146000 • e-mail: info@adler-healthcare.com • www.adlerhealthcare.com